

Client: Peptora, LLC
help@peptorapeptides.com
855-811-0112

Sample received: **03/30/26**
Analysis conducted: **04/03/26**

Compound:	Ipamorelin	CAS:	170851-70-4
Batch/Lot #:	INIPA10031226	Formula:	C38H49N9O5
Appearance:	White lyophilized powder	Mol Wt:	711.9 g/mol


Method: Qualitative and Quantitative chemical analysis by Ultra High Performance Liquid Chromatography with Mass Spectrometry

Pubchem CID: 9831659

[Ipamorelin](#) | [C38H49N9O5](#) | [CID 9831659 - PubChem](#)


	Specification	Result	
Identity Test:	Ipamorelin	Ipamorelin	Conforms
Quantity:	10 mg	10.62 mg	Conforms
Purity:	>99%	99.41%	Conforms

LC-MS Chromatogram and Retention Time Analysis




Full Scan Mass Spectrometry Analysis

PP-Ipamorelin-10mg-01 #700 RT: 3.05 AV: 1 NL: 2.50E8
T: FTMS + p ESI Full ms [300.0000-1300.0000]



Calibration Curve



Analysis Performed by
Dr. Victor Ugur
Analytical Chemist
contact@bioregen.com

COA #12718
Security Key **PEPTORA01**
bioregen.com/verify



Victor Ugur

Limulus Amebocyte Lysate (LAL) – Report of Analysis

Test Performed	Endosafe - PTS LAL Assay
Test Date	04/08/2026
Analyst	VICTOR
Instrument	Endosafe® nexgen-PTS™ (SN: 25201216)

Sample Information

Company Name	Peptora, LLC
Sample Tested	Ipamorelin
Sample Lot Number	INIPA10031226
Sample Concentration	10mg
Vials/Units Received	1

Reagent Information

Cartridge Type	Endotoxin- LAL (Charles River)
Cartridge Lot #	2602741388
Cartridge Calibration Code	518334799106
Detection Range	0.05 EU/mL – 5.0 EU/mL
Spike Value	1.28 EU/mL
Spike Recovery	116% (Pass)
Test Suitability	Pass

Endotoxin Quantitation Results

Sample ID	Replicates	Measured EU/mL	Specification	Pass/Fail
Ipamorelin	1, 2	<0.05 EU/mL	≤ 0.50 EU/mL	Pass

Interpretation:

Endotoxin level in **Ipamorelin** was measured using chromogenic LAL assay with detection at 405 nm. The result of **<0.05 EU/mL** indicates that the endotoxin content is **below** the acceptable threshold of 0.5 EU/mL. Sample **PASSES** the endotoxin limit test.



Analyst Signature: _____

Date: 04/08/2026

To verify the authenticity of this report, please contact **BioRegen** at contact@bioregen.com or reference the **Report ID** and **Validation Key** provided above.

Alternatively, scan the QR code below to access the official record on our website.

